


# Lab automation built for the future

With modular design and intelligent software, the possibilities for innovation are limitless.




#### Introduction

# You Innovate, We Automate

Automating your scientific workflow doesn't have to be intimidating or take a huge investment of time and resources. With the Opentrons Flex, you have an entire team of scientific and engineering experts backing you up. With modular hardware, unmatched usability, comprehensive support, and available in the configuration that meets your needs. Run your experiments your way, just more of them....a lot more

We're making automation accessible. How will you Flex it?



# Opentrons FLEX Solution and and a series of the series s

#### Modular system evolves with your needs

Increase throughput, change pipettes or add modules to automate more complex workflows

#### Maximize walkaway time and throughput

Run high-throughput workflows up to 96 samples

#### Accessible to all

Easy to use, no-code tools and extensive support

#### High performance, priced for innovation

Without compromising quality

## Flex Fits Your Needs



#### **On-Instrument Guided Setup**

Dynamic deck maps match your protocol and guide you through protocol setup

#### Protocol Selection on the Flex

Send protocols from the desktop app to Flex for easy browsing and initiation

#### **Automated Calibration**

One-touch automated pipette and gripper positional calibration

#### **Protocol Development Your Way**

No-code Protocol Designer app, Protocol Library, Python API, or design services to build your protocol

#### **Third-Party Integrations**

Open source platform enables easy integrations

#### **Flex Specifications**

#### NUMBER OF DECK SLOTS

12 main deck slots, plus 4 additional slots for staging tips and labware

#### **VOLUME RANGE**

 $1-1000 \mu L$ 

#### **ACCURACY\***

- +/- 8.00% systematic error at  $1 \mu L$
- +/- 7.00% random error at 1 µL

#### PIPETTE CONFIGURATION

Any two Flex 1-Channel or 8-Channel Pipettes, or either Flex 96-Channel Pipette alone

#### LABWARE COMPATIBILITY

Any automation-compatible, ANSI/SLAS standard footprint well plates, plus a wide range of other labware, including reservoirs and tubes up to 50 mL conical

#### CONNECTIVITY

USB, Ethernet, and Wi-Fi

#### **DIMENSIONS**

87 cm x 69 cm x 84 cm (W, D, H)

\*Measured with the Opentrons Flex 1-Channel Pipette (1–50  $\mu L)$  using a 50  $\mu L$  tip

# Flex Fits Your Application



Opentrons Flex workstations come with the Flex robot and accessories, pipettes and gripper, on-deck modules, and labware you need to automate your chosen application.\* You can configure each workstation with either Flex 1-and 8-Channel Pipettes or the Flex 96-Channel Pipette to suit your scale. The workstations are reagent- and kit-agnostic to minimize changes to your current workflow.



#### ♦ Automate up to 96 samples

#### **Choose from these** workstations covering some of our most popular applications:

#### OPENTRONS FLEX HIGH-THROUGHPUT **NGS WORKSTATION**

Automation of library preparation, enabling full walkaway processing for up to 96 samples

#### **OPENTRONS FLEX PROTEOMICS** WORKSTATION

Automation of proteomics workflows including ELISA, mass spec sample prep, and more!

#### OPENTRONS FLEX NGS WORKSTATION

For end-to-end NGS library prep

#### OPENTRONS FLEX PCR WORKSTATION

For PCR setup and thermocycling

#### **OPENTRONS FLEX NUCLEIC ACID EXTRACTION WORKSTATION**

For nucleic acid extraction

#### **OPENTRONS FLEX PROTEIN PURIFICATION WORKSTATION**

For small-scale protein purification using magnetic beads

# Flex Fits Your Expectations



| OPENTRONS FLEX<br>1-CHANNEL PIPETTES | OPENTRONS FLEX<br>8-CHANNEL PIPETTES | OPENTRONS FLEX<br>96-CHANNEL PIPETTE |
|--------------------------------------|--------------------------------------|--------------------------------------|
| 1–50 μL                              | 1–50 μL                              | 1–200 μL                             |
| 5–1000 μL                            | 5–1000 μL                            | 5–1000 μL                            |

Flex pipettes use air displacement technology to offer highly accurate pipetting from 1 to 1000  $\mu$ L. Smart sensors support automatic calibration, real-time positioning, and error detection. Choose the pipettes that fit your workflow — the gantry supports any two Flex 1-Channel or 8-Channel Pipettes, or the Flex 96-Channel Pipette alone — and swap them out yourself when your workflow changes.

|   | VOLUME (μL) | SPECIFICATION |        |
|---|-------------|---------------|--------|
|   |             | D (%)         | CV (%) |
| OPENTRONS FLEX 1-CHANNEL PIPETTE (1–50 μL)*†        | 1           | 8.00          | 7.00   |
|   | 10          | 1.50          | 0.50   |
|   | 50          | 1.25          | 0.40   |
| OPENTRONS FLEX<br>1-CHANNEL PIPETTE<br>(5-1000 μL)‡ | 5           | 5.00          | 2.50   |
|   | 50          | 0.50          | 0.30   |
|   | 200         | 0.50          | 0.15   |
|   | 1000        | 0.50          | 0.15   |

<sup>\*</sup>Flex 8-Channel and 96-Channel Pipette specifications available from <a href="mailto:info@opentrons.com">info@opentrons.com</a>

 $<sup>^{\</sup>dagger}All$  volumes tested with 50  $\mu L$  tip

 $<sup>^{\</sup>dagger}5$  and 50  $\mu$ L volumes tested with a 50  $\mu$ L tip; 200  $\mu$ L volume tested with a 200  $\mu$ L tip; 1000  $\mu$ L volume tested with a 1000  $\mu$ L tip

## Flex Fits Your Workflow

Opentrons hardware modules are designed, tested and built by Opentrons. They work seamlessly with Opentrons robots and software. Every Opentrons Flex Module that sits on the deck comes with a caddy to secure it into the deck, putting labware close to the deck surface and keeping cables out of the way.



#### **GRIPPER**

The Opentrons Flex Gripper attaches to the gantry and is capable of moving labware across the deck and on or off modules such as the Magnetic Block or Thermocycler. It can also access the additional slots to the side of the main deck to replenish tips or labware.



#### THERMOCYCLER MODULE

The Opentrons Thermocycler Module GEN2 is a fully automated on-deck thermocycler, providing hands-free PCR in a 96-well plate format. The heated lid and disposable seal fit tightly over the plate, ensuring efficient sample heating and minimal evaporation.



#### TEMPERATURE MODULE

The Opentrons Temperature Module GEN2 is an easy-to-use hot and cold plate that can maintain constant temperatures from 4 to 95 °C. Each module comes with aluminum blocks to hold labware, including 1.5 mL tubes, 2 mL tubes, PCR plates and strips, and most flat bottom well plates.



#### MAGNETIC BLOCK

The Opentrons Flex Magnetic Block is a 96-well magnetic well plate holder containing static ring magnets. Use it with the Flex Gripper and skirted 96-well PCR or deep-well plates to facilitate quick and efficient separation using magnetic beads.



#### **HEATER-SHAKER MODULE**

The Opentrons Heater-Shaker Module provides on-deck heating and orbital shaking of microplates. The module can be heated to 95  $^{\circ}$ C, and can shake samples from 200 to 3000 rpm. It is compatible with a number of deep-well and 96-well plates.

#### Flex Modules



#### **STACKER**

The Opentrons Flex Stacker is a modular tip and plate storage and delivery system designed to enhance the throughput of your laboratory workflows while conserving valuable deck space. Each stacker can accommodate up to 7 tip racks, 48 PCR plates, or 16 deep well plates. Integrate up to four stackers without significantly increasing the instrument's footprint.



#### ABSORBANCE PLATE READER

The Opentrons Flex Absorbance Plate Reader is a fully automated solution for absorbance-based assays. Capable of detecting absorbance in the 450–1000 nm range, it is optimized for a broad array of applications, including protein quantification, sample normalization, cell viability assays, and bacterial growth monitoring.



#### HEPA/UV

Keep your work contamination free with the Opentrons Flex HEPA/UV Module. With UV-C sterilization and filters exceeding HEPA standards, rest assured your samples are safe from microorganisms, nucleic acids, and nucleases.



#### **DECK EXPANSION**

The Opentrons Flex Deck Expansion Set allows you to utilize up to 4 more positions on the Opentrons Flex worktable. The expansions slots are installed to the right of the worktable, acting as positions A4, B4, C4, and D4 and accepting plates, tube racks, and other labware that support SBS dimensions.



#### **WASTE CHUTE**

The Gripper moves and drops tip racks and labware into the Waste Chute and a trash receptacle below to increase walkaway capacity. The Waste Chute accepts plates, tip racks and tips maximizing your throughput.



# We're making automation accessible. How will you Flex it?

#### Need more information?

- For general sales information, email <a href="mailto:info@opentrons.com">info@opentrons.com</a>
- Want to speak to a sales rep? Mon Fri, Global: +1 (646) 212-3703
   UK/EU: +44 7588 611072

Trademarks: Opentrons Flex®, Opentrons drop logo, (Opentrons Labworks, Inc.). Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

APRIL 2025 © OPENTRONS 2025, ALL RIGHTS RESERVED