




Opentrons Flex[®] Stacker

5x the throughput without 5x the space




How will you Flex it?

Maximize your throughput with the Opentrons Flex Stacker

The Opentrons Flex Stacker is a modular tip and plate storage and delivery system designed to enhance the throughput of your laboratory workflows while conserving valuable deck space. Easily installable onto the Opentrons Flex expansion slots, each stacker can accommodate up to 7 tip racks, 48 PCR plates, or 16 deep well plates. With the capability to integrate up to four stackers without significantly increasing the instrument's footprint, the Flex Stacker enables efficient processing of up to 96 NGS samples in a single run, making it an ideal solution for high-throughput genomics and proteomics applications



Genomics




Proteomics



High-throughput
screening

Increased Deck Space

The Opentrons Flex Stacker allows you to fit up to five times more labware into your workflow without increasing the instrument's footprint, enabling high-throughput processing while preserving valuable deck space.



Seamless and Modular Installation

Designed for seamless integration, the Flex Stacker can be easily installed onto the Opentrons Flex expansion slots, allowing you to add up to four stackers as your throughput needs grow.

Effortless Programming

With Opentrons' intuitive software solutions, incorporating the Flex Stacker into your protocols is straightforward, enabling efficient automation of complex workflows.



The Opentrons Flex Stacker in Action!

To demonstrate the effectiveness of the Opentrons Flex Stacker, we successfully automated two popular DNA library prep chemistries, that use mechanical shearing and enzymatic fragmentation respectively, at 96-sample throughput.

Key Findings

The Opentrons Flex Stacker supports walkaway automation of popular DNA library preparation kits, resulting in high-yielding, consistent libraries.




Figure 1. Deck layout for automating library preparation for 96 samples using IDT xGEN MC¹ and Twist EF 2.0² kits on the Opentrons Flex High-Throughput NGS Workstation

Row A: (A1) On Deck Thermocycler Module with Sample Plate, (A2-A3) Tip Racks. (A4) Flex Stacker Module with Tip Racks

Row B: (B2) Waste Reservoir, (B3) Reagent Plate 2, (B4) Flex Stacker Module with Tip Racks

Row C: (C1) Temperature Module with Reagent Plate 1, (C2) Reagent Reservoir, (C3) Plate Lids, (C4) Deck Extension with Final Sample Plate

Row D: (D1) Tip Rack, (D2) Cleanup Plate 1, (D3) Waste Chute, (D4) Deck Extension with Cleanup Plate 2

Table 1. Sequencing Results depicting the GC Content, Percent Passing Filter, Quality >30 and Percent Mapped for Lambda Samples prepared with IDT xGen MC and Human Genomic Samples prepared with Twist EF 2.0, on a 2 x 150 MiSeq Run.

IDT xGen MC			Twist EF 2.0	
	Lambda	CV	Human Genomic	CV
GC content	49.21%	2.27%	41.85%	0.12%
% PF	98.55%	0.47%	99.40%	0.00%
Q30	94.44%	0.56%	94.74%	0.12%
% Mapped	98.73%	1.02%	100.00%	0.00%




Figure 2. IDT xGen MC Library Preparation (A) Sample Plate Layout, (B) Library Fragment Sizes, 2% E-Gel, E-Gel Sizing DNA Ladder, (C) Library Yield, Flex and Manual Comparison.




Figure 3. Twist Library Preparation Enzymatic Fragmentation Kit 2.0 (A) Sample Plate Layout, Row G and H Tested in a Checkerboard pattern, (B) Library Fragment Sizes, Library Fragment Sizes 2% E-Gel, E-Gel Sizing DNA Ladder, representative Positive and Negative Samples, (C) Library Yield, Positive and Negative Control Comparison.

Specifications

Installation Capacity	4
Compatible Labware	Opentrons Flex Tip Racks, PCR plates, flat bottom plates, deep well plates
Labware Capacity	7 tip racks, 48 PCR / flat bottom plates, 16 deep well plates
Material	Aluminum / Polycarbonate
Dimensions (mm)	957.5 x 106 x 194.5
Weight	18kg / 40lbs
Electrical Input	6.4A 115VAC, 3.2A 230VAC
Electrical Output	48V, 2.7A



Get hands-on support, services and training with our Flex Service offerings.

- **12-Month Factory Warranty**
Included with all robot purchases
- **Installation**
Robot unboxing and setup
- **Onsite or Remote Custom Protocol Development**
Consult your sales rep

Need more information?

- For general sales information, email info@opentrons.com
- For services and installation, email support@opentrons.com
- Want to speak to a sales rep? Mon – Fri, Global: +1 (646) 212-3703
UK/EU: +44 758 861 1072