




Technical Data Sheet

Opentrons Flex™ Tips

Opentrons Flex 1000 µL, Filter Tips - SKU: 991-00106



Tip volume	1000 µL
Tip length	95.50 mm
Packaging quantity	96 pcs/box
Material	Polypropylene
Temperature Range	<ul style="list-style-type: none">• Melting Point: 160° C• Heat Deflection: 70-100° C
Color	Clear
Sterility	<ul style="list-style-type: none">• Assurance Level: SAL 10-6• Free of human DNA, PCR inhibitors, RNase/DNase, pyrogens, and endotoxins

