





# Technical Data Sheet

## Opentrons Flex™ Tips



### Opentrons Flex 1000 µL Tips - SKU: 991-00103



<b>Tip volume</b>	1000 µL
<b>Tip length</b>	95.50 mm
<b>Packaging quantity</b>	96 pcs/box
<b>Material</b>	Polypropylene
<b>Temperature Range</b>	<ul style="list-style-type: none"><li>• Melting Point: 160° C</li><li>• Heat Deflection: 70-100° C</li></ul>
<b>Color</b>	Clear
<b>Sterility</b>	<ul style="list-style-type: none"><li>• Assurance Level: SAL 10-6</li><li>• Free of human DNA, PCR inhibitors, RNase/DNase, pyrogens, and endotoxins</li></ul>

