


Technical Data Sheet

Pipette Accuracy



Flex 50 µl Single Channel Pipette	Tip Size	Volume Range	µL	Spec	
				% d	% CV
	50 µl	1-50	1	8%	7%
	50 µl	1-50	10	1.50%	0.50%
	50 µl	1-50	50	1.25%	0.4%



Flex 1000 µl Single Channel Pipette	Tip Size	Volume Range	µL	Spec	
				% d	% CV
	50 µl	5-50	5	5%	2.5%
	50 µl	5-50	50	0.5%	0.30%
	200 µl	5-200	200	0.5%	0.15%
	1000 µl	5-1000	1000	0.5%	0.15%



Flex 50 µl Multi Channel Pipette	Tip Size	Volume Range	µL	Spec	
				% d	% CV
	50 µl	1-50	1	10%	8%
	50 µl	1-50	10	2.5%	1%
	50 µl	1-50	50	1.25%	0.6%



Flex 1000 µl Multi Channel Pipette	Tip Size	Volume Range	µL	Spec	
				% d	% CV
	50 µl	5-50	5	8%	4%
	50 µl	5-50	50	2.5%	0.6%
	200 µl	5-200	200	1%	0.25%
	1000 µl	5-1000	1000	0.7%	0.15%



Flex 96 Channel Pipette	Tip Size	Volume Range	µL	Spec	
				% d	% CV
	50 µl	5-50	5	10%	5%
	50 µl	5-50	50	2.5%	1.25%
	200 µl	5-200	200	1.5%	1.25%
	1000 µl	5-1000	1000	1.5%	1.5%