



# Opentrons Flex MagBead Protein Purification Workstation

For scaling up and fully automating magnetic bead-based protein purification workflows

For automation of 8 to 96 samples  
1- & 8-Channel Pipettes

\$45,750 ~~\$50,830~~

SAVE 10%

Part number: 991-00361

For automation of up to 96+ samples  
96-Channel Pipette

\$52,000 ~~\$57,780~~

SAVE 10%

Part number: 991-00362

- Robot and Pipettes
- Tips and Labware
- Magnetic Block
- Heater-Shaker and Caddy
- Waste Chute



With the Flex Protein Purification Workstation, you can automate your small-scale protein purification and proteomics sample processing at the scale you need to increase efficiency, reduce errors, and save hands-on time. This Opentrons Protein Purification Workstation uses magnetic beads for purification of proteins.

This workstation can be used for the following applications:

- Enzymatic and Chemical Lysis
- SP3 Sample Preparation for LCMS
- BCA and Bradford Assays
- Protein Labeling
- SDS and Native PAGE Sample Preparation

ITEM	1- & 8-CHANNEL WORKSTATION	96-CHANNEL WORKSTATION
Opentrons Flex Robot	1	1
Opentrons Flex Gripper GEN1	1	1
Opentrons Flex 1-Channel (5-1000 µL)	1	
Opentrons Flex 8-Channel (5-1000 µL)	1	
Opentrons Flex 96-Channel (5-1000 µL)		1
Opentrons Flex Magnetic Block GEN1	1	1
Heater-Shaker Module with Opentrons Flex Caddy	1	1
Opentrons Flex Waste Chute GEN2	1	1
Opentrons Flex Deck Expansion Set (4 count)		1
Opentrons Flex Filter Tips, 200 µL	20 racks (20 refills) Equivalent to 3,840 tips	20 racks (20 refills) Equivalent to 3,840 tips
Opentrons Flex Filter Tips, 1000 µL	20 racks (20 refills) Equivalent to 3,840 tips	20 racks (20 refills) Equivalent to 3,840 tips
4-in-1 Tube Rack Set	1	1
NEST 2.0 mL Microcentrifuge Tubes (500 count)	1	1
NEST 15 mL Centrifuge Tube (500 count)	1	1
NEST 2.0 mL Deep Well Plate, V Bottom (50 count)	1	1
Opentrons Tough PCR Plate, Full Skirt (25 count)	1	1
NEST 12-well Reservoirs, 15 mL (50 count)	1	1
NEST 1-well Reservoirs, 290 mL (50 count)	1	1
Remote Protocol Development	1	1

- 12-Month Factory Warranty**  
 Included with all robot purchases
- Opentrons Care**  
 Remote service contract with 1, 2 or 3 year options
- Opentrons Care Plus**  
 On-site service contract with 1, 2 or 3 year options
- Installation**  
 Consult your sales rep
- Upgrade Options Available**  
 Consult your sales rep



**Customize your Flex Workstation to fit your workflow by adding on any of the following modules at 10% savings.**

- Opentrons Temperature Module
- Cerillo Plate Reader
- Opentrons Flex HEPA/UV Module

**NEED MORE INFORMATION?**

- For general sales information, email [info@opentrons.com](mailto:info@opentrons.com)
- For services and installation, email [support@opentrons.com](mailto:support@opentrons.com)
- Want to speak to a sales rep? Call us between 9:00 AM – 6:00 PM EST, Mon – Fri, +1 (908) 869-8907