



Opentrons Flex Nucleic Acid Extraction Workstation

For scaling up and fully automating magnetic bead-based nucleic acid extraction

For automation of up to 48 samples
1- & 8-Channel Pipettes

\$44,300 ~~\$49,222~~

SAVE 10%

Part number: 991-00183


For automation of up to 96 samples
96-Channel Pipette

\$52,700 ~~\$58,555~~

SAVE 10%

Part number: 991-00184

- Robot and Pipettes
- Heater-Shaker Module and Magnetic Block
- Tips and Labware



With the Flex Nucleic Acid Extraction Workstation, you can automate your NAE workflow at the scale you need to increase efficiency, reduce errors, and save hands-on time. Flex enables you to automate your nucleic acid isolation and purification workflows using many leading magnetic bead-based reagent systems on the market.

This workstation can be used for the following applications:

- Sample lysis via heat and/or agitation
- Magnetic bead mixing and separation

ITEM	1- & 8-CHANNEL WORKSTATION	96-CHANNEL WORKSTATION
Opentrons Flex Robot	1	1
Opentrons Flex Gripper GEN1	1	1
Opentrons Flex 1-Channel (5-1000 µL)	1	
Opentrons Flex 8-Channel (5-1000 µL)	1	
Opentrons Flex 96-Channel (5-1000 µL)		1
Opentrons Flex Magnetic Block GEN1	1	1
Heater-Shaker GEN1 with Opentrons Flex Caddy	1	1
Heater-Shaker PCR Plate Adapter	1	1
Heater-Shaker Deep Well Plate Adapter	1	1
Opentrons Flex Waste Chute GEN1	1	1
Opentrons Flex Deck Expansion Set		1
Opentrons Flex Filter Tips, 200 µL	80 racks (20 refills) Equivalent to 19,200 tips	80 racks (20 refills) Equivalent to 19,200 tips
Opentrons Flex Filter Tips, 1000 µL	160 racks (20 refills) Equivalent to 38,400 tips	280 racks (20 refills) Equivalent to 57,600 tips
NEST 12-Well Reservoirs, 15 mL (50 count)	1	1
NEST 1-Well Reservoirs, 290 mL (50 count)	1	1
Opentrons Tough 0.2 mL 96-Well PCR Plate, Full Skirt (25 count)	2	2
NEST 2 mL Deep Well Plate, V Bottom (50 count)	1	2
4-in-1 Tube Rack Set	1	1
Remote Protocol Development	1	1

Get hands-on support, services and training with our Flex Service offerings.

- **12-Month Factory Warranty**
Included with all robot purchases
- **Opentrons Protect**
12-month extended warranty
- **Opentrons Care**
Remote service contract with 1, 2 or 3 year options
- **Opentrons Care Plus**
On-site service contract with 1, 2 or 3 year options
- **Installation**
Consult your sales rep
- **Onsite or Remote Custom Protocol Development**
Consult your sales rep

NEED MORE INFORMATION?

- For general sales information, email info@opentrons.com
- For services and installation, email support@opentrons.com
- Want to speak to a sales rep? Call us between 9:00 AM – 6:00 PM EST, Mon – Fri, +1 (908) 869-8907



Customize your Flex Nucleic Acid Extraction Workstation to fit your workflow by adding on any of the following modules at 10% savings.

- Opentrons Temperature Module GEN2 with Opentrons Flex Caddy
- Opentrons Flex Deck Expansion Set (4 count) (1- & 8-Channel Workstation Only)
- Opentrons Flex HEPA/UV Module