



Opentrons Flex PCR Workstation

For scaling up and fully automating PCR setup and thermocycling

For automation of 8 to 96 samples
1- & 8-Channel Pipettes

\$42,750 ~~\$47,500~~

SAVE 10%

Part number: 991-00356


For automation of 96+ samples*
96-Channel Pipette

\$48,000 ~~\$53,333~~

SAVE 10%

Part number: 991-00357

- Robot, Pipettes, and Waste Chute
- Temperature module and optional Thermocycler
- Tips and Labware



With the Flex PCR Workstations, you can automate your PCR setup and thermocycling at the scale you need to increase efficiency, reduce errors, and save hands-on time. Configure your workstation with our on-deck Thermocycler for end-to-end automation, or with no Thermocycler for automated PCR setup with off-deck thermocycling.

This workstation can be used for the following applications:

- PCR setup
- Thermocycling

*Workstation can be used to automate 96 samples of PCR setup and on-deck thermocycling, or 192 samples of PCR setup (i.e., with off-deck thermocycling).

ITEM	1- & 8-CHANNEL WORKSTATION	96-CHANNEL WORKSTATION
Opentrons Flex Robot	1	1
Opentrons Flex Gripper GEN1	1	1
Opentrons Flex 1-Channel (1-50 µL)	1	
Opentrons Flex 8-Channel (1-50 µL)	1	
Opentrons Flex 96-Channel (5-1000 µL)		1
Temperature Module GEN2 with Opentrons Flex Caddy	1	1
Opentrons Flex Waste Chute GEN2	1	1
Opentrons Flex Filter Tips, 50 µL	20 racks (20 refills) Equivalent to 3,840 tips	20 racks (20 refills) Equivalent to 3,840 tips
NEST 2.0 mL Microcentrifuge Tubes (500 count)	1	1
NEST 12-Well Reservoirs, 15 mL (50 count)	1	1
NEST 1-Well Reservoirs, 290 mL (50 count)		1
Opentrons Tough 0.2 mL 96-Well PCR Plate, Full Skirt (25 count)	1	1
4-in-1 Tube Rack Set	1	1
Remote Protocol Development	1	1

Get hands-on support, services and training with our Flex Service offerings.

- 12-Month Factory Warranty**
 Included with all robot purchases
- Opentrons Care**
 Remote service contract with 1, 2 or 3 year options
- Opentrons Care Plus**
 On-site service contract with 1, 2 or 3 year options
- Installation**
 Consult your sales rep
- Onsite or Remote Custom Protocol Development**
 Consult your sales rep



Customize your Flex PCR Workstation to fit your workflow by adding on any of the following modules at 10% savings.

- Thermocycler GEN2 with Opentrons Flex Caddy
- 2nd Temperature Module GEN2 with Opentrons Flex Caddy
- Opentrons Flex Deck Expansion Set (4 count)
- Opentrons Flex HEPA/UV Module

NEED MORE INFORMATION?

- For general sales information, email info@opentrons.com
- For services and installation, email support@opentrons.com
- Want to speak to a sales rep? Call us between 9:00 AM – 6:00 PM EST, Mon – Fri, [+1 \(908\) 869-8907](tel:+19088698907)