



Opentrons Flex Plasmid Prep Workstation

For scaling up and fully automating magnetic bead-based plasmid extraction and purification

For automation of 8 to 96 samples
1- & 8-Channel Pipettes

\$44,500 \$49,415

SAVE 10%

Part number: 991-00187

- Robot and Pipettes
- Tips and Labware
- Magnetic Block
- Heater-Shaker
- Waste Chute
- Remote Protocol Development



The Flex Plasmid Prep Workstation is built to support magnetic bead based plasmid extraction and purification workflows. This workstation is equipped with high volume pipettes, heater-shaker module, and magnetic block to accommodate for most bead based chemistry. Protocol development is Included in each workstation to ensure the Flex is ready to run soon after installation.

This workstation can be used for the following applications:

- Bacterial cell lysis
- Magnetic bead separation
- Plasmid washing and elusion

ITEM	1- & 8- CHANNEL WORKSTATION
Opentrons Flex Robot	1
Opentrons Flex 1-Channel (5-1000 µL)	1
Opentrons Flex 8-Channel (5-1000 µL)	1
Flex Gripper	1
Flex Magnetic Block	1
Flex Waste Chute	1
Flex Heater-Shaker with Flex Caddy	1
Opentrons Flex Filter Tips, 200 µL	20 racks (20 refills) Equivalent to 3,840 tips
Opentrons Flex Filter Tips, 1000 µL	20 racks (20 refills) Equivalent to 3,840 tips
NEST 12-Well Reservoirs, 15 mL (50 count)	1
NEST 1-Well Reservoirs, 290 mL (50 count)	1
NEST 2.0 mL Microcentrifuge Tubes (500 count)	1
NEST 2.0 mL 96-Well Deep Well Plate, V Bottom (50 count)	1
Opentrons Tough 0.2 mL 96-Well PCR Plate, Full Skirt (25 count)	1
4-in-1 Tube Rack Set	1
Remote Protocol Development	1

Get hands-on support, services and training with our Flex Service offerings.

- 12-Month Factory Warranty**
 Included with all robot purchases
- Opentrons Protect**
 12-month extended warranty
- Opentrons Care**
 Remote service contract with 1, 2 or 3 year options
- Opentrons Care Plus**
 On-site service contract with 1, 2 or 3 year options
- Installation**
 Consult your sales rep
- Upgrade Options Available**
 Consult your sales rep

NEED MORE INFORMATION?

- For general sales information, email info@opentrons.com
- For services and installation, email support@opentrons.com
- Want to speak to a sales rep? Call us between 9:00 AM – 6:00 PM EST, Mon – Fri, +1 (908) 869-8907



Customize your Flex Plasmid Prep Workstation to fit your workflow by adding on any of the following modules at 10% savings.

- Temperature Module and Caddy
- Deck Expansion Set