



Opentrons Flex Prep Workstation

For scaling up and fully automating general liquid handling such as sample transfers and reagent aliquoting.

For automation of 8 to 96 samples
1- & 8-Channel Pipettes

\$26,400 ~~\$29,235~~

SAVE 10%

Part number: 991-00187

For automation of 96+ samples*
96-Channel Pipette

\$32,200 ~~\$35,778~~

SAVE 10%

Part number: 991-00188

- Robot and Pipettes
- Tips and Labware
- Remote Protocol Development



The Flex Prep Workstation is a simple and affordable starter bundle suitable for simple pipetting workflows such as sample transfers/duplication, reagent aliquoting and plate stamping. This customizable workstation offers optional add-on modules or pipettes to support flexible application needs. Protocol development is included in each workstation to ensure the Flex is ready to run soon after installation.

This workstation can be used for the following applications:

- Sample and Reagent Aliquoting
- Serial Dilutions
- Cherry Picking
- Normalization and Pooling

ITEM	1- & 8-CHANNEL WORKSTATION	96-CHANNEL WORKSTATION
Opentrons Flex Robot	1	1
Opentrons Flex 1-Channel (5-1000 µL)	1	
Opentrons Flex 8-Channel (5-1000 µL)	1	
Opentrons Flex 96-Channel (5-1000 µL)		1
Opentrons Flex Filter Tips, 200 µL	20 racks (20 refills) Equivalent to 3,840 tips	20 racks (20 refills) Equivalent to 3,840 tips
Opentrons Flex Filter Tips, 1000 µL	20 racks (20 refills) Equivalent to 3,840 tips	20 racks (20 refills) Equivalent to 3,840 tips
NEST 12-Well Reservoirs, 15 mL (50 count)	1	1
NEST 1-Well Reservoirs, 290 mL (50 count)	1	1
NEST 2.0 mL Microcentrifuge Tubes (500 Count)	1	1
4-in-1 Tube Rack Set	1	1
Remote Protocol Development	1	1

Get hands-on support, services and training with our Flex Service offerings.

- 12-Month Factory Warranty**
 Included with all robot purchases
- Opentrons Protect**
 12-month extended warranty
- Opentrons Care**
 Remote service contract with 1, 2 or 3 year options
- Opentrons Care Plus**
 On-site service contract with 1, 2 or 3 year options
- Installation**
 Consult your sales rep
- Upgrade Options Available**
 Consult your sales rep



Customize your Flex Prep Workstation to fit your workflow by adding on any of the following modules at 10% savings.

- Opentrons Flex Gripper
- Opentrons Flex Heater-Shaker + Caddy
- Opentrons Flex Waste Chute
- Opentrons Flex HEPA/UV Module

NEED MORE INFORMATION?

- For general sales information, email info@opentrons.com
- For services and installation, email support@opentrons.com
- Want to speak to a sales rep? Call us between 9:00 AM – 6:00 PM EST, Mon – Fri, [+1 \(908\) 869-8907](tel:+19088698907)